


## SpikeLog128 Specification



This document provides specifications and instructions that are specific to the SpikeLog128 neural logger. General instructions for use of Deuteron's neural loggers and related software is available on Deuteron's website [www.deuterontech.com](http://www.deuterontech.com)

### General description

SpikeLog128 is a single board neural logger that provides 128 channels of neural recording. It also logs data from a "9-axis" motion sensor. SpikeLog128B hardware can also connect to external sensors and log data coming from them.

The main blocks of this logger are:

- Four 32-channel digitizing amplifier array ICs
- A digital radio IC that uses sub GHz ISM band radio frequencies
- A MicroSD memory card; up to 512HGz is presently supported
- A microcontroller
- A logic block designed to move bulk data rapidly between digital sections of the logger
- Two battery chargers
- A 9-axis motion sensor IC

The neural inputs connect via two 70-pin connectors. Electrode assemblies if use with this logger must either be terminated with a matching pair of 70 pin loggers, or, more usually, the logger plugs into an adapter that connects between the electrode assembly and a matching pair of 70-pin connectors.

#### Basic capabilities

Function	Properties	Detail
<b>Presentation</b>	Single circuit board, vertical orientation	No standard housing
<b>Neural recording</b>	64 channels	32000 samples per second
<b>Motion sensor</b>	Optional. "9-Axis"; 3D-accelerometer, 3D-gyroscope and 3D-magnetometer	Gyro and accelerometer: 1000 samples per second per axis. Magnetometer: 110 samples per second per axis
<b>Data recording</b>	To removable MicroSD card	Cards up to 512GB are supported
<b>Wired data streaming</b>	None	
<b>Radio data streaming</b>	Preview only	Limited to about 100Kbytes per second

## Specifications

### File formats

If configured for neural data only, Deuteron's "Flat" file format can be used. Systems that support motion sensor recording, audio recording, or both use deuteron's "block" file format.

### Neural Recording

Function	Properties	Detail
Signal range	10mV p-p	
Digital resolution	0.2 $\mu$ V	
Random noise	2.4 $\mu$ V	For 7kHz bandwidth
Input capacitance	14pF	
Preamplifier bandwidth	Low limit :0.2 to 500Hz High limit: 200Hz to 10kHz	Fully software selectable
Analog filters	Low-pass: 3 <sup>rd</sup> order. High-pass: 1 <sup>st</sup> order	
Sampling rate	32000 samples per second, each channel	Fixed sampling rate
Connector	Pair of HiroseDF40C-70DP-0.4V(51)	Spacing between connector center lines is 25.00mm
Reference channels	Ground is the default reference, but one of seven channels, or ground, in each bank can be selected to be a reference channel.	There are four banks of 32 signals. The voltage of the selected reference channel is subtracted from all the channels in that bank before the signals are amplified

## Motion Sensor


Note: Motion sensor functions are available only in software versions that support block file format.

Function	Properties	Detail
<b><u>Accelerometers</u></b>		
Number of axes	3	
Data rate for each axis	1000 ± 10	Samples per second
Measurement ranges	±2g, ±4g, ±8g, ±16g	4 software-selectable ranges
Resolution	61μg, 123μg, 246μg, 492μg	16-bit, according to range
Noise PSD	300μg / sqrt Hz	
<b><u>Gyroscopes</u></b>		
Number of axes	3	
Data rate for each axis	1000 ± 10	Samples per second
Measurement ranges	250°/s, 500°/s, 1000°/s, 2000°/s	4 software-selectable ranges
Resolution	0.076°/s, 0.015°/s, 0.031°/s, 0.061°/s	16-bit, according to range
Noise PSD	0.01 (°/s) / sqrt Hz	
<b><u>Magnetometers</u></b>		
Number of axes	3	
Data rate for each axis	111 ± 3	Samples per second
Measurement range	±4.8	mT
Resolution	0.6μT	14 bit
Uncorrected zero error	300	μT


## Physical and electrical

Function	Properties	Detail
Dimensions	30 x 30 x 4mm	Except for Micro-SD connector
Mass	4.4g	Including microSD card
Battery voltage	3.7V	Designed for any lithium polymer protected cell above 100mAh
Charger current	300mA	Suitable for cells 300-1200mAh, other currents available on request
Charger connector	Micro USB	5V supply required
SD card socket	Push-in Push out	
Current consumption	102mA	Typical when recording, motion sensor off.
Memory card capacities	8GB – 512GB	

Co-ordinate axes for Accelerometer and Gyroscope




Co-ordinate axes for Accelerometer and Gyroscope



## Rear View

For front view see first page of this specification



## Available adapters

### SpikeLog128 for Cereport Adapter

Cereport is an elastomer-based mount developed in academia and sold commercially by Blackrock. Below is a view of the Cereport adapter assembled to the SpikeLog128 board, as well as two views of the complete logger for use on monkeys enclosed in a strong plastic housing that holds the logger, adapter, battery and radio antenna. The memory card can be inserted and removed via a slot in the housing. Full documentation for the Cereport adapter is provided in a separate document.

