

## SpikeLog32 Specifications



This document provides specifications and instructions that are specific to the SpikeLog32 neural logger. General instructions for use of Deuteron's neural loggers and related software is available on Deuteron's website [www.deuterontech.com](http://www.deuterontech.com)

### General description

SpikeLog32 is a single board neural logger that provides 32 channels of neural recording. It also can record audio and ultrasonic signals using its on-board microphone or from an external transducer. A 9-axis motion sensor can be used to record motion and direction.

## Basic capabilities

Function	Properties	Detail
<b>Presentation</b>	Single circuit board, vertical orientation	No standard housing
<b>Neural recording</b>	32 channels	32000 samples per second
<b>Audio / ultrasonic recording</b>	Optional. 1 channel, up to 80kHz bandwidth	50K, 100K or 200k samples per second
<b>Motion sensor</b>	Optional. "9-Axis"; 3D-accelerometer, 3D-gyroscope and 3D-magnetometer	Gyro and accelerometer: 1100 samples per second per axis. Magnetometer: 110 samples per second per axis
<b>Data recording</b>	To removable MicroSD card	Cards up to 512GB are supported
<b>Wired data streaming</b>	Optional: USB streaming of all data, file download, impedance measurement	See USB interface documentation
<b>Radio data streaming</b>	Preview only	Limited to about 100Kbytes per second

## Specifications

### File formats

If configured for neural data only, Deuteron's "Flat" file format can be used. Systems that support motion sensor recording, audio recording, or both use deuteron's "block" file format.

### Neural Recording

Function	Properties	Detail
<b>Signal range</b>	10mV p-p	
<b>Digital resolution</b>	0.2 $\mu$ V	
<b>Random noise</b>	2.4 $\mu$ V	For 7kHz bandwidth
<b>Input capacitance</b>	15pF	
<b>Preamplifier bandwidth</b>	Low limit :0.2 to 500Hz High limit: 200Hz to 10kHz	Fully software selectable
<b>Analog filters</b>	Low-pass: 3 <sup>rd</sup> order. High-pass: 1 <sup>st</sup> order	
<b>Sampling rate</b>	32000 samples per second, each channel	Fixed sampling rate
<b>Connector</b>	Omnetics A79025-001	Four corner pins are ground pins
<b>Reference channels</b>	Ground is the default reference, but six of the 32 input pins channel can optionally be selected as a reference channel.	The voltage of the selected reference channel is subtracted from all of the channels before the signals are amplified

## Physical and electrical

Function	Properties	Detail
<b>Dimensions</b>	23 x 16 x 3mm	
<b>Mass</b>	2.3g	Including microSD card
<b>Battery voltage</b>	3.7V	Designed for any lithium polymer protected cell above 40mAh
<b>Current consumption</b>	45mAh	Typical when recording, motion sensor off.
<b>Memory card capacities</b>	8GB – 512GB	

## Motion Sensor

Note: Motion sensor functions are available only in software versions that support block file format.


Function	Properties	Detail
<b><u>Accelerometers</u></b>		
Number of axes	3	
Data rate for each axis	1000 ± 10	Samples per second
Measurement ranges	±2g, ±4g, ±8g, ±16g	4 software-selectable ranges
Resolution	61µg, 123µg, 246µg, 492µg	16-bit, according to range
Noise PSD	300µg / sqrt Hz	
<b><u>Gyroscopes</u></b>		
Number of axes	3	
Data rate for each axis	1000 ± 10	Samples per second
Measurement ranges	250°/s, 500°/s, 1000°/s, 2000°/s	4 software-selectable ranges
Resolution	0.076°/s, 0.015°/s, 0.031°/s, 0.061°/s	16-bit, according to range
Noise PSD	0.01 (°/s) / sqrt Hz	
<b><u>Magnetometers</u></b>		
Number of axes	3	
Data rate for each axis	111 ± 3	Samples per second
Measurement range	±4.8	mT
Resolution	0.6µT	14 bit
Uncorrected zero error	300	µT

## Audio Recording

Function	Properties	Detail
<b>On-board sensor</b>	Knowles SPU0410LR5H-QB	Sensitive beyond 80kHz
<b>Audio bandwidth. Upper limit</b>	80kHz	
<b>Audio bandwidth, Lower limit</b>	300Hz	
<b>Audio sampling rates</b>	50kHz, 100kHz, 200kHz	
<b>Audio sensitivity</b>	High sensitivity: 60 $\mu$ Pa/bit (5dBA per bit level)	16 bit range is 80dBA

Spectral response of microphone in audio and ultrasonic ranges as characterized by Knowles.

**Typical Free Field Response  
Normalized to 1kHz**



**Preliminary Ultrasonic Free Field  
Response Normalized to 1kHz**

